

Supplementary Table 2. Regents used in experiments.

Regents	Clone	Application	Supplier
Alexa Fluor® 647 Mouse Anti-Human TIM-3	565558	Flowcytometry	BD Pharmingen
PE Mouse Anti-Human CD279	557946	Flowcytometry	BD Pharmingen
FITC anti-human CD69	310904	Flowcytometry	BioLegend
Live/Dead Fixable Red Dead Stain	L10120	Flowcytometry	Thermofisher Scientific
APC anti-human CD8	344722	Flowcytometry	BioLegend
CD8 MicroBeads, human	130-045-201	Cells sorting	Miltenyi
Dynabeads Human T-Activator CD3/CD28	11131D	T cell activation	Thermofisher Scientific
Alexa Fluor® 488 anti-human Ki-67	558616	Flowcytometry	BD Pharmingen
APC anti-human IL-2	500310	Flowcytometry	BioLegend
Alexa Fluor® 488 anti-human IFN- γ	502515	Flowcytometry	BioLegend
PE anti-human CD8 ,	344706	Flowcytometry	Biolegend
Alexa Fluor® 647 Mouse IgG1 κ Isotype Control	565571	Flowcytometry	BD Pharmingen
PE Mouse IgG1, κ Isotype Control	555749	Flowcytometry	BD Pharmingen
FITC Mouse IgG1 κ Isotype Ctrl Antibody	488108	Flowcytometry	Biolegend
Alexa Fluor 488 Mouse IgG1 κ Isotype Ctrl	400129	Flowcytometry	Biolegend
APC Rat IgG2a, κ Isotype Ctrl	400512	Flowcytometry	Biolegend
Rabbit anti-CD63	10606D	Immunofluorescence	Thermofisher Scientific